


ADPT-EIB-36-QC-HS-36
Input: QuickClip 72
Output: Omnetics 44 pin X 2



ADPT-EIB-16-QC-HS-16
Input: QuickClip 16
Output: Omnetics 26 pin



ADPT-HS-18-PUSH-PIN
Input: 2x5 833 Mill-Max
Output: Omnetics 26 pin




ADPT-HS-18-Versa
Input: Mill-Max 18 pin
Output: Omnetics 26 pin



ADPT-EIB-54-HS-72-QC
Input: Mill-Max 14/16/18 socket
Output: QuickClip 72



ADPT-EIB-54-Kopf-HS-36
Input: Mill-Max 14/16/18 socket
Output: Omnetics 44 pin X 2



ADPT-HS-36-GMR
Input: Omnetics 30 socket
Output: Omnetics 44 pin

