





ADPT-PED-128-HS-36

Input: Blackrock Pedestal
Output: Omnetics 44 X 4




ADPT-HS-36-QC-N2T-32

Input: Omnetics 36 socket
Output: QuickClip 36



ADPT-HS-16-QC-N2T-16

Input: Omnetics 18 socket
Output: QuickClip 16



ADPT-HS36-N2T-64-Inline

Input: SamTex 4x10 X 2
Output: Omnetics 44 X 2



ADPT-HS36-N2T-32A

Input: SamTec 4x10
Output: Omnetics 44



ADPT-HS36-N2T-32

Input: SamTec 4x10
Output: Omnetics 44



ADPT-HS18-N2T-16R

Input: Omnetics 18 socket
Output: Omnetics 28 pin




ADPT-HS27-N2T-16

Input: 16 pin DIP
Output: Mill-Max 36



ADPT-HS18-N2T-16

Input: Omnetics 18 socket
Output: Omnetics 28 pin



ADPT-HS18-N2T-16A

Input: 16 pin DIP
Output: Omnetics 28 pin