


FreeLynx™

FreeLynx is **the** wireless acquisition system for freely moving subjects. With 32 to 256 *configurable* channels of wide bandwidth wireless neural recording and real-time signal processing with reliable transmission up to 10 meters, depending on environment.



- 30 kHz, 16 bit data acquisition
- Dual 9-axis accelerometer
- Omni-directional video tracking LEDs
- Setup options for weight distribution on subject
- Easy swap batteries with minimal recording downtime and fast reboot
- Switch between onboard microSD or wireless transmission at any time

Not intended for human use.