

# LabLynx™

A full-featured Neuralynx data acquisition system for wide bandwidth neural recording, perfect for new or established labs. With wireless transmission to the acquisition computer and internal battery power option, the system can be truly portable and easily moved between experiment locations. The LabLynx connects up to 4 digital multiplexing headstages for low equipment costs and weight.



- 16 - 256 channels, 16 bit resolution, up to 30 kHz sample rate
- 10 meters transmission range, depending on environment
- Rechargeable battery or AC power source
- Digital I/O for experiment control
- 0.1 Hz - 8 kHz frequency range from each electrode